



SHORTCUT SHEET

SAY-IT! SAM TABLET



Modify Page:

- Open to page you want to change
- Press book icon at top (menu)
- Click on Setup (tools at bottom)
- Click on page setup
- Click on edit current page
- Click on edit page cells
- Choose button to change
- Click on Edit Cell
- Click in box to put in label (should be short). Type in label and then click on close.
- Click on box to put in message (can be as long as needed). Type in message and then click on close.
- (If keyboard does not pop up for either label or message, then click on the other box and it should pop up)
- Click on empty square button on right side to put a picture in
- Search—if it can't find it, type in a new word to search for.
- Choose picture
- Close—sometimes it will bring you back to that page and you will need to re-choose the picture.
- Check message line to make sure extra letters/words didn't appear!
- On Action pull down menu, click on clear/add message/speak
- Done, Yes to save

Reset Say it Sam:

Press reset, calendar, envelope, and power until you get black screen, then let go. It will go through a series of restarts that takes about 5 minutes. When you get to the loading sample page, go to menu, setup, layout, load new, and click on your layout.

http://www.words-plus.com/website/pdf_files/SAM2CommManual.<http://www.words-plus.com/website/pdf_files/SAM2CommManual.pdf>pdf

If you ever have other issues that you cannot fix, contact www.words-plus.com or call them at 800-869-8521