

SHORTCUT SHEET

TECH/TALK



The Tech/Talk is an augmentative communication device that records message and plays them back using keys with graphic overlays.

TO USE

1. Turn it on (turn volume knob on the back clockwise, it clicks)
2. Check the **Level** knob to select the Level you want to use.
3. Press the squares to hear the message.

RECORDING

1. Messages have to be short - 1 -2 seconds each
2. For best speech quality, turn volume down low, record with a strong voice then playback with the volume up higher.
3. Make sure you are prepared and know exactly what you want each location to say (write it down before you start recording)
4. Select the **Level** you want to change.
5. Press and hold the **Record** button on the back while at the same time hold down the square on which you are recording the message.
6. About 6 – 8 inches away from the machine speak the message loud and clearly.
7. Release both buttons immediately after you have finished speaking (so you don't record dead sound)
8. The microphone is located on the back panel or on the front underside depending on how old your machine is.
9. **Test** your new message by pressing the message square you just recorded on.

TROUBLESHOOTING

If the device stops working it is probably because the batteries have run out. It uses 4 AA batteries. You may need a small Phillips screwdriver to open the battery compartment (or strong fingernails) depending on the model you have.

PRODUCT INFORMATION

Contact Information: AMDi
Advanced Multimedia Devices, Inc
200 Frank Road
Hicksville, NY 11801
(888) 353-2634

Pricing: \$395 to \$645
Price varies based on the number of levels (6, 8, or 12).

Edited (Created) by Lori O'Hare, SWAAAC, 4/19/2004

For more SWAAAC Shortcut Sheets, go to www.swaacac.org, click on **AT Loan Bank** and then click on **Shortcut Sheets**. Would you like a Shortcut Sheet that doesn't exist? Contact the WestAC office at (970) 248-0876 or carrie.foster@uchsc.edu.